


The *NANOStation II*

High-power scanning probe microscopy

The *NANOStation II* combines optical microscopy and scanning probe microscopy (SPM) in a single, optimized set-up. The combination of the tried and trusted S.I.S. *ULTRAObjective* scanning probe microscope and a powerful optical microscope allows unequalled productivity during the high-resolution inspection of surfaces. The pre-selection of interesting structures for the SPM is greatly simplified through the use of a high-quality optical microscope. The base of the *NANOStation II* is a very sturdy microscope stand which has been optimized against vibrations and thermal drift. Equipped with a high-power Zeiss optical microscope and the S.I.S. *ULTRAObjective* scanning probe microscope, the *NANOStation II* is an exceptionally versatile inspection microscope.

SURFACE IMAGING SYSTEMS




Measuring modes

All *ULTRAObjective* measuring modes are available

Top: Optical image of an etched GaAs crystal (DIC)
3D image of the aforementioned etching, measured in non-contact mode with the *NANOStation II*

Far left: *NANOStation II* with special table to inspect CDs

Left: *NANOStation II* with special holder to inspect contact lenses in a physiological solvent



Applications



The *NANOStation II* can be easily adapted to special requirements and customers' wishes.

Example 1, CD and DVD inspection: A special table was developed, which uses a suction component to lock the sensitive discs in place without stress.

Example 2, Contact Lens Inspection: A holder securing the lenses and permitting measurements at defined points was developed to allow inspection of contact lenses in a liquid.

Far left: 2D image of a collagen fibre, measured in non-contact mode with the *NANOStation II*

Left: Cross-section through the image to the left



Specifications of the *NANOStation II*-SPM-System

Scan range:	20 μm x 20 μm x 3 μm 40 μm x 40 μm x 4 μm 80 μm x 80 μm x 5 μm 100 μm x 100 μm x 10 μm 200 μm x 200 μm x 10 μm hardware linearized scan motion in X-Y-direction (optional in Z-direction)	Input channels:	max. 4 simultaneous
Noise level:	0,1 nm rms in vertical direction (Z)	External inputs:	max. 3 high speed with 16 bit resolution
Lateral accuracy:	typically within 1%, closed loop scanning	Image size:	freely selectable, from 128 to 1024 pixels, even rectangular sizes
Scan speed:	typ. 1 to 10 Hz	Processing:	internal 32 bit DSP, typ. 50 MHz
Detection principle:	fiber optical interferometry, noise level < 0,01 nm rms	Computer interface:	USB (standard universal serial bus)
Tips:	silicon tips, various types	Operating system:	MS-Windows 2000®
Tip change:	adjustment free	Microscope:	Zeiss Axiotech optional with bright/dark field or differential interference contrast (DIC)
Digital input resolution:	16 bit A/D	Positioning:	manual translation stage 25 mm x 25 mm other sizes available on request
Digital output resolution:	16 bit D/A	Weight:	approximate 50 kg
Output voltage:	\pm 165 V, with 2 μV rms	Material:	granite